

## Beta Spectrum

**Goal:** to investigate the spectrum of  $\beta$  rays emitted by a  $^{137}\text{Cs}$  source. The instrument used is a so-called  $180^\circ$  magnetic spectrometer that separates  $\beta$  rays of different momenta by the fact that they follow different paths in a magnetic field.

### 1. Introduction

The continuous spectrum of  $\beta$  rays emitted from a radioactive source has been explained by Fermi<sup>1</sup> [FER34]. Introductions into the theory of  $\beta$  decay can be found in many texts<sup>2, 3, 4, 5</sup>. Aside from the Fermi decay, there exists also a process called 'Internal Conversion' (of nuclear excitation energy) which leads to *discrete* lines in the  $\beta$  spectrum (ref. 2, p.396, ref. 5, p.362). A lot of interesting information on internal conversion is contained in a PhD thesis<sup>6</sup> of a former IU student. Here, it is assumed that you have learned enough about  $\beta$  decay to be able to answer the following questions: What determines the shape of the continuous spectrum? What is a Kurie plot? Can you explain the Internal Conversion process?

In this experiment we study the  $\beta$  decay of the  $^{137}\text{Cs}$  nucleus that includes both, the Fermi decay and the internal conversion. You can find a lot of detailed information on this decay in ref. 7. Try to answer the following: What is the half-life of  $^{137}\text{Cs}$ ? What is the highest energy  $\beta$  rays that are emitted from  $^{137}\text{Cs}$ ? What role does the nucleus  $^{137}\text{Ba}$  play in all of this? What are the expected energies of the conversion lines? Draw a decay scheme with the relevant details using information from ref. 7.

### 2. Equipment

The principle of operation of magnetic  $\beta$  spectrometers is explained in ref.2, p. 52-63. Some important relations between the parameters of the instrument and the quantities to be measured are also given there.


In the 1950's the measurement of  $\beta$  activity in nuclei was an important part of the nuclear research in the IU Physics Department, and a group headed by L.M. Langer (†2000) made many important contributions to the field. The spectrometer used in our lab is a piece of history of physics at IU. It is a somewhat smaller version of the original device<sup>8, 9</sup>, built in 1948. Cook's thesis<sup>9</sup> contains a detailed description of the apparatus and many experimental details, and you should definitely read it. The smaller, newer device is also described in the literature<sup>10</sup>. A copy of ref. 10 is enclosed.

In the magnetic field of the spectrometer  $\beta$  rays from a source are guided on a semi-circular path to a slit in front of a detector. The pole pieces are cleverly shaped, providing a radial field dependence such that about 8% of the  $\beta$  rays emitted by the source are focused onto the exit slit. At a given magnetic field  $B$  (in Gauss),  $\beta$  rays of a certain momentum  $p$  (in keV/c) are


transported around the semi-circular path with a mean bending radius  $\rho$  (in cm). In our instrument,  $\rho=15\text{cm}$ . The three quantities are related by the Lorentz force:

$$p = 0.3 \cdot B \cdot \rho \tag{1}$$

The range of accepted momenta is defined by a slit of variable opening  $s$  (see fig.1). The slit width governs the momentum resolution  $\Delta p$  of the spectrometer. Using a typical slit width of  $s = 1 \text{ mm}$ , you should estimate  $\Delta p$  as the difference in momentum between two  $\beta$  particles with bending radii  $\rho$  and  $\rho+s/2$ . In measuring a spectrum, it is clearly important that the magnetic field strength is known well<sup>11</sup>. Instead of using eq.1, it is more accurate to calibrate the relation between momentum and field by observing  $\beta$ 's of a known momentum, e.g., from the K conversion line in  $^{137}\text{Ba}$  (see ref.12).



The detector consists of a piece of plastic scintillator (fig.1) glued to a light guide that is coupled to a RCA 8575 photomultiplier. Pulses from the photomultiplier are sent to a discriminator which selects only pulses above an adjustable threshold for counting. This is



**Figure 2.** Energy loss by electrons in air, Al, and Pb. To suppress the large variation in  $dE/dx$  arising from the number of electrons of the material, the quantity  $\rho^{-1}(dE/dx)$  is plotted. Solid lines are for collisions; dashed lines are for radiation. For additional tabulated data on energy losses, see L. Pages et al., *Atomic Data* 4, 1 (1972).

needed in order to discriminate against the small, but frequent noise pulses. Since the pulse height also depends on the energy of the  $\beta$  particles, it is important to set the threshold carefully.

The inside of the spectrometer is evacuated and a thermocouple gauge measures the pressure in the spectrometer chamber. The question arises: how good does the vacuum have to be? In order to answer this, you should calculate the energy loss along the semi-circular path of  $\beta$ 's with momenta between 0.1 and 1.0 MeV/c and compare this figure to the momentum resolution  $\Delta p$ . The energy loss for  $\beta$  particles in various materials is given in fig.2.

### 3. Measurements

- the performance of photomultipliers is affected by magnetic fields: does the fringe field from the spectrometer have an effect on your measurement?
- how to set the threshold of the discriminator?
- measure the magnetic field  $B$  with the Hall probe and the flip coil.

- the magnet current  $I$  can be varied (useful range is from 20mA to 100mA): establish the connection between  $I$  and  $B$ ; are there hysteresis effects?.
- calibrate the magnetic field using the known rigidity of the K conversion line ref. 2, p.227.
- measure the continuous part of the spectrum; construct a Kurie plot; determine the endpoint energy.
- how many conversion lines are visible? determine their energies.
- determine the conversion coefficient  $\kappa$  of the K line.
- determine the intensity ratios between conversion lines.

## 4. References

- <sup>1</sup> E. Fermi, *Z.f.Physik* **88** (1934) 161
- <sup>2</sup> K. Siegbahn, *Beta- and Gamma-Ray Spectroscopy*, North-Holland, Amsterdam 1955. (QC771.S57)
- <sup>3</sup> W.E. Burcham, *Elements of Nuclear Physics*, Longman, London 1979, p.295ff. (QC776.B85)
- <sup>4</sup> H. Frauenfelder and E.M. Henley, *Subatomic Physics*, Prentice-Hall, Englewood 1974, p.273ff. (QC776.F845)
- <sup>5</sup> E. Segré, *Nuclei and Particles*, Benjamin, Reading 1977, p.410ff. (QC776.S4)
- <sup>6</sup> G.A. Graves, PhD thesis, Physics Department, Indiana University, 1952 (QC1000.G776)
- <sup>7</sup> C.M. Lederer and V.S. Shirley, *Table of Isotopes*, Wiley, New York 1978.
- <sup>8</sup> L.M. Langer and C.S. Cook, *Rev. Sci. Instr.* **19**, 257 (1948)
- <sup>9</sup> C.S. Cook, PhD thesis, Physics Department, Indiana University, 1948 (QC1000.C77)
- <sup>10</sup> J.A. Bruner and F.R. Scott, *A High-Resolution Beta-Ray Spectrometer*, *Rev. Sci. Instr.* **21**, 545 (1950).
- <sup>11</sup> L.M. Langer and F.R. Scott, *Rev. Sci. Instr.* **21**, 522 (1950)
- <sup>12</sup> L.M. Langer and R.D. Moffat, *Phys. Rev.* **78**, 74 (1950)