



Univerza v Ljubljani



SiPMs for Čerenkov imaging


Peter Križan

University of Ljubljana and J. Stefan Institute


Trends in Photon Detectors in Particle Physics and Calorimetry,
Trieste, June 2-4, 2008





Contents



Photon detectors for Ring Imaging Cherenkov counters
Example: proximity focusing RICH for the Belle upgrade
SiPMs as single photon counters (uniformity, timing)
Detection of Cherenkov photons with SiPMs
SiPMs for PET
Summary




Measuring Cherenkov angle



Proximity focusing RICH


Idea: transform the **direction** into a **coordinate** →
ring on the detection plane
→ **Ring Imaging Cherenkov**



RICH with a focusing mirror




Photon detection in RICH counters



RICH counter: measure photon impact point on the photon detector surface

→ detection of **single** photons with

- sufficient **spatial resolution**
- **high efficiency** and **good signal-to-noise ratio**
- over a **large area** (square meters)





Special requirements:

- **Operation in magnetic field**
- **High rate capability**
- **Very high spatial resolution**
- **Excellent timing (time-of-arrival information)**



Belle upgrade – side view




Two new particle ID devices, both RICHes:

Barrel: **Time Of Propagation (TOP)** counter or **focusing DIRC**


Endcap: **proximity focusing RICH**



Endcap: Proximity focusing RICH



K/ π separation at 4 GeV/c:
 $\theta_c(\pi) \sim 308$ mrad ($n = 1.05$)
 $\theta_c(\pi) - \theta_c(K) \sim 23$ mrad



For single photons: $\delta\theta_c(\text{meas.}) = \sigma_0 \sim 14$ mrad,
typical value for a 20mm thick radiator and 6mm PMT pad size

Per track:
$$\sigma_{\text{track}} = \frac{\sigma_0}{\sqrt{N_{pe}}}$$

Separation: $[\theta_c(\pi) - \theta_c(K)] / \sigma_{\text{track}}$


$\rightarrow 5\sigma$ separation with $N_{pe} \sim 10$




Beam tests

pion beam (π^2) at KEK

π




Photon detector: array of 16 H8500 ('flat pannel') PMTs



Aerogel

Clear rings, little background



June 4, 2008

TPDPPC




Photon detectors for the aerogel RICH



Needs:

- Operation in high magnetic field (1.5T)
- High efficiency at $\lambda > 350\text{nm}$
- Pad size $\sim 5\text{-}6\text{mm}$






Candidates:

- MCP PMT (Burle 85011)
- large area HAPD of the proximity focusing type (R+D)
- SiPM?



Photon detector candidate: Large active area HAPD




R&D project in collaboration with Hamamatsu.

Long development time, now working test samples.

First beam test results →

→NIMA, RICH07 proc.


TPDPPC08, Trieste







Photon detector candidate: BURLE/Photonis MCP-PMT

BURLE 85011 microchannel plate (MCP) PMT: multi-anode PMT with two MCP steps




- good performance in beam and bench tests, NIMA567 (2006) 124
- very fast →
- R+D: ageing









Photon detector candidate: BURLE/Photonis MCP-PMT



BURLE 85011 microchannel plate (MCP) PMT: time resolution after time walk correction




Tails can be significantly reduced by:

- decreased photocathode-MCP distance and
- increased voltage difference

→ NIMA, RICH07 proc.



SiPM as photon detector?



Can we use SiPM (Geiger mode APD) as the photon detector in a RICH counter?

- +immune to magnetic field
 - +high photon detection efficiency, single photon sensitivity
 - +easy to handle (thin, can be mounted on a PCB)
 - +potentially cheap (not yet...) silicon technology
 - +no high voltage
-
- very high dark count rate (100kHz – 1MHz) with single photon pulse height
 - radiation hardness





Can such a detector work?



Experience from HERA-B RICH:
successfully operated in a high
occupancy environment (up to
10%).


Need **>20** photons per ring (had
 ~ 30) for a reliable PID.

HERA-B RICH event







Can such a detector work?




MC simulation of the counter response: assume 1mm^2 active area SiPMs with 0.8 MHz (1.6 MHz, 3.2 MHz) dark count rate, 10ns time window

K identification efficiency at 1% π missid. probability



For different background levels






SiPMs as photon detectors?

SiPM is an array of APDs operating in Geiger mode. Characteristics:


- low operation voltage $\sim 10\text{-}100\text{ V}$
- gain $\sim 10^6$
- peak PDE up to 65%(@400nm)

$$\text{PDE} = \text{QE} \times \epsilon_{\text{geiger}} \times \epsilon_{\text{geo}}$$


- ϵ_{geo} – dead space between the cells
- time resolution $\sim 100\text{ ps}$
- works in high magnetic field
- dark counts $\sim \text{few } 100\text{ kHz/mm}^2$
- radiation damage (p,n)



Hamamatsu MPPC: S10362-11






Surface sensitivity for **single** photons




- 2d scan in the focal plane of the laser beam ($\sigma \approx 5 \mu\text{m}$)
- intensity: on average $\ll 1$ photon \rightarrow single photons
- Selection: single pixel pulse height, in 10 ns TDC window

S137 (CPTA/Photonique)








Surface sensitivity for single photons 2




E407 (Pulsar/MEPHI)






Surface sensitivity for single photons 3





Hamamatsu MPPCs

H050C




H025C








Surface sensitivity for single photons 4



Hamamatsu MPPCs

H100C









Time resolution: time walk correction





<< 1 photon

uncorrected TDC






corrected TDC






Time resolution: blue vs red




	E407	S137	H100C	H050C	H025C
σ_{red} (ps)	127	182	145	212	154
σ_{blue} (ps)	97	151	136	358	135

• $\sigma \approx 100$ ps

• $\sigma_{\text{red}} > \sigma_{\text{blue}}$



Expected number of photons for aerogel RICH




with multianode PMTs or SiPMs(100U), and
aerogel radiator: thickness 2.5 cm, $n = 1.045$
and transmission length (@400nm) 4 cm.



$$N_{\text{SiPM}}/N_{\text{PMT}} \sim 5$$

Assuming 100% detector
active area

Never before tested in a RICH
where we have to detect single
photons. ← Dark counts have
single photon pulse heights
(rate 0.1-1 MHz)



SIPMs: Cosmic test setup



- 6 Hamamatsu SiPMs (=MPPC) of type 100U (10x10 pixels with 100 μ m pitch), background \sim 400kHz
- signals amplified (ORTEC FTA820),
- discriminated (EG&G CF8000) and
- read by multihit TDC (CAEN V673A) with 1 ns / channel


June 4, 2008

TPDPPC0

To be published in NIM A in \sim 3 weeks



SiPM: Cherenkov angle distributions for 1ns time windows







-6 ns

-5 ns

-4 ns

-3 ns



+2 ns

+3 ns

+4 ns





+5 ns

Cherenkov photons appear in the expected time windows →

First Cherenkov photons observed with SiPMs!



SiPM Cherenkov angle distribution




Fit function is a combination of

- a background (quadratic) and
- a signal (Gaussian).


Only scale parameters are free

→ SiPMs give 5 x more photons than PMTs per photon detector area – in agreement with expectations

To be published in NIM A in ~3 weeks





SIPMs: improving signal/noise




Improve the signal to noise ratio:

- Reduce the noise by a narrow (few ns) time window
- Increase the number of signal hits per single sensor by using light collectors and by adjusting the pad size to the ring thickness

Light collector with reflective walls





or combine a lens and mirror walls





Light collection: improve signal to noise ratio

Machined from a plastic plate
(HERA-B RICH lens material).



Light guide should be as close as possible to the SiPM surface
(now: air gap + epoxy layer)









Cherenkov photons with light collectors




First attempt: use the top of a blue LED






- ★ Yield increase in agreement with the expectations
- ★ Further improvements possible by
 - Using SiPMs with a reduced epoxy protective layer
 - using a better light collector





Detector module design







Light guide geometry optimisation





Light Guide Acceptance / (d and out)






Light guide geometry optimisation



d (mm)	out (mm)	accept. (%)
3.0	1.48	51.6
3.1	1.45	54.0
3.2	1.41	55.7
3.3	1.38	57.8
3.4	1.34	59.2
3.5	1.31	61.0
3.6	1.27	62.6
3.7	1.24	63.1
3.8	1.20	64.4
3.9	1.16	64.4
4.0	1.13	64.9
4.1	1.09	64.3
4.2	1.06	63.8
4.3	1.02	62.8
4.4	0.99	61.8
4.5	0.95	60.5
4.6	0.92	58.5
4.7	0.88	56.4
4.8	0.85	54.6
4.9	0.81	51.9


SiPM = 0.8, M = 3.3, d = 5.0 | gap(y,z) = (0.0, 0.0) | $\theta = 30.0$ Thu May 8 14:02:15 2008



Detector module – final version for beam tests this week at KEK

SiPMs: array of 8x8 SMD mount
Hamamatsu S10362-11-100P
with 0.3mm protective layer

Light guides



2cm


SiPMs + light guides


June 4, 2008

Photon detector for the beam test


64 SiPMs

20mm





Open question: sensitivity to neutrons




However: SiPMs are sensitive to neutron irradiation (dark count rate starts increasing drastically after $\sim 10^{10}$ neutrons/cm²)

- We have to measure the neutron flux in the relevant detector region: **calibrated Si** diodes mounted in the spectrometer since January, leave it there for ~ 1 year, extract and determine the integrated flux (fluence)
- We will also mount a few **SiPMs** in the proper place in the spectrometer, register their performance during running.



Summary



Single Cherenkov photons were observed for the first time with SiPM in a RICH counter using cosmic rays

Small size light guides were designed, machined and attached to the SiPMs.

SiPMs are a promising candidate for photon detection in future RICH counters


Plan:

Further study of different light collection systems


Test a larger array

Open issue: influence of neutrons on the counter performance

Explore other applications of the device → PET



We also work on a PET module...




Test a PET module with:


4x4 array of LYSO crystals ($4.5 \times 4.5 \times 20(30) \text{ mm}^3$)

SiPMs: Photonique $2.1 \times 2.1 \text{ mm}^2$


SiPMs: Photonique
 $2.1 \times 2.1 \text{ mm}^2$



LYSO crystals
 $4.5 \text{ mm} \times 4.5 \text{ mm}$





We also work on a PET module 2




Some tests with Na22 in coincidence with a 4x4
LYSO+MAPMT module

Best: $\sim 9\%$ (rms) energy resolution

Shown: one of the early results.






Next step...




Test a PET module with:



4x4 array of LYSO crystals (4.5 x 4.5 x 20(30) mm³)

16 SiPMs (Photonique 2.1x2.1 mm²)





Next step...




4x4 SiPMs

Reading out a
4x4 array of
LYSO crystals

Module under test...




Back-up slides







Planar entry window






Spherical entry window



Spherical entry window, reflective sides



Efficiency vs. angle of incidence α





Light guide	d/a	R/a	$\alpha_{\min}, \alpha_{\max}$	$I(-60^\circ, 60^\circ)$
Planar entry	3.4	-	$-24^\circ, 24^\circ$	64%
Sph. entry	1.6	2.0	$-35^\circ, 35^\circ$	66%
Reflective sides	2.4	2.6	$-44^\circ, 44^\circ$	69%



TOF capability of a RICH

With a fast photon detector (MCP PMT), a proximity focusing RICH counter can be used also as a **time-of-flight counter**.


Time difference between π and K \rightarrow



For time of flight: use Cherenkov photons emitted in the **PMT window**




TOF capability: window photons



Expected number of detected Cherenkov photons emitted in the PMT window (2mm) is **~15**

→ Expected resolution **~35 ps**



TOF test with pions and protons at 2 GeV/c.


Distance between start counter and MCP-PMT is 65cm

→ In the real detector ~2m

→ 3x better separation




Time-of-flight with photons from the PMT window



Benefits: Čerenkov threshold in glass (or quartz) is much lower than in aerogel.

Aerogel: kaons (protons) have **no** signal below 1.6 GeV (3.1 GeV): identification in the **veto** mode.



Threshold in the **window:** π K p



Window: threshold for kaons (protons) is at ~ 0.5 GeV (~ 0.9 GeV): \rightarrow **positive identification** possible.



Photon detector candidate: MCP-PMT

BURLE 85011 MCP-PMT:

- multi-anode PMT with two MCP steps
- 25 μm pores
- bialkali photocathode
- gain $\sim 0.6 \times 10^6$
- collection efficiency $\sim 60\%$
- box dimensions $\sim 71\text{mm}$ square
- 64(8x8) anode pads
- pitch $\sim 6.45\text{mm}$, gap $\sim 0.5\text{mm}$
- active area fraction $\sim 52\%$



- Tested in combination with multi-anode PMTs


- $\sigma_g \sim 13 \text{ mrad}$ (single cluster)
- number of clusters per track $N \sim 4.5$
- $\sigma_g \sim 6 \text{ mrad}$ (per track)
- $\rightarrow \sim 4 \sigma \pi/K$ separation at 4 GeV/c

- 10 μm pores required for 1.5T
- collection eff. and active area fraction should be improved
- aging study should be carried out


Light collection: single vs multi channel



Multichannel device+imaging light collection system: Has a very limited angular acceptance




Single channel: combine a lens and mirror walls






Beam test: Cherenkov angle resolution and number of photons



NIM A521(2004)367; NIM A553(2005)58

Beam test results with 2cm thick aerogel tiles:


>4 σ K/ π separation



→ Number of photons has to be increased.




Radiator with multiple refractive indices



How to increase the number of photons without degrading the resolution?

normal


→ stack two tiles with different refractive indices: “focusing” configuration






→ focusing radiator



Focusing configuration – data



4cm aerogel single index



2+2cm aerogel

